



## Product Specification

**Blood Gas Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)**

**Page 1**



### Product description

|                     |   |
|---------------------|---|
| Order number        | 05.1146   |
| Product description | Blood Gas Monovette® , preparation: calcium-balanced lithium heparin, nominal volume: 1 ml, cap white/orange, label/print: transparent/orange, connection: Luer (m), (LxØ) without cap: 66 x 11 mm, 50 piece(s)/case, sterile, 50 piece(s)/case |

### Product characteristics

|                           |                                  |
|---------------------------|----------------------------------|
| Type of collection        | Arterial and venous              |
| Type of preparation       | calcium-balanced lithium heparin |
| Preparation concentration | 50 IU/ml blood                   |
| Label/ Print              | with plastic label               |
| Colour of print/label     | transparent/orange               |
| Graduation                | fill mark                        |
| Application               | blood gas analysis               |
| Connection                | Luer (m)                         |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Blood Gas Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)**

**Page 2**

### Size

|                                  |        |
|----------------------------------|--------|
| Sample volume                    | 1 ml   |
| Diameter                         | 11 mm  |
| Length including cap             | 82 mm  |
| Length excluding cap             | 66 mm  |
| Length including cap and plunger | 100 mm |

### Material & colours

|                         |                                   |
|-------------------------|-----------------------------------|
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | white/orange                      |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

### Purity & certification

|                  |                                |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified     | CE                             |
| Purity standard  | sterile                        |
| Sterilisation    | Electron irradiation           |
| Batched          | yes                            |

### Packaging

|                               |        |
|-------------------------------|--------|
| Minimum order qty.            | 500    |
| Type of smallest subpackaging | case   |
| Piece(s) / inner box          | 50     |
| Piece (s) / outer case        | 500    |
| Piece(s) / pallet             | 56000  |
| Depth of box                  | 145 mm |
| Width of box                  | 108 mm |
| Height of box                 | 63 mm  |
| Depth of case                 | 333 mm |
| Width of case                 | 303 mm |
| Height of case                | 125 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Blood Gas Monovette® , calcium-balanced lithium heparin, 1 ml, cap white/orange, connection: Luer (m)**

**Page 3**

|                   |               |
|-------------------|---------------|
| Case volume       | 0.0126 cbm    |
| Weight of product | 0.0043 kg     |
| Weight of case    | 2.69 kg       |
| EAN of inner box  | 4038917247221 |
| EAN case          | 4038917247214 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-22

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**