



Product Specification
Ventilator for Blood Gas Monovette® ,
transparent/orange, length: 24.5 mm



Product description

Order number	14.1148
Product description	Ventilator for Blood Gas Monovette® , transparent/orange, length: 24.5 mm, accessories for the Blood Gas Monovette® , 500 piece(s)/case

Product characteristics

Application	accessories for the Blood Gas Monovette®
Example of use	Ventilation of the Blood Gas -Monovette®
Suitable for	Blood Gas Monovette®

Size

Length of product	24.5 mm
-------------------	---------

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent/orange

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Ventilator for Blood Gas Monovette® ,
transparent/orange, length: 24.5 mm

Page 2

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	168000
Depth of box	191 mm
Width of box	129 mm
Height of box	129 mm
Depth of case	288 mm
Width of case	218 mm
Height of case	150 mm
Case volume	0.0094 cbm
Weight of product	0.0009 kg
Weight of case	1.34 kg
EAN of inner box	4038917085465
EAN case	4038917075312

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2023-02-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT