



**Product Specification**  
**Blood gas capillaries, calcium-balanced lithium  
heparin, 100 µl**

**Page 1**



**Product description**

Order number	19.930.100
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 100 µl, (LxØ): 100 x 2.05 mm, 250 piece(s)/container

**Product characteristics**

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	80 IU/ml blood
Application	blood gas analysis

**Size**

Sample volume	100 µl
Diameter	2.05 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Blood gas capillaries, calcium-balanced lithium  
heparin, 100 µl**

**Page 2**

**Material & colours**

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

**Packaging**

Minimum order qty.	250
Type of smallest subpackaging	container
Piece(s) / inner box	250
Piece (s) / outer case	16000
Piece(s) / pallet	960000
Depth of case	398 mm
Width of case	358 mm
Height of case	175 mm
Case volume	0.0249 cbm
Weight of product	0.0003 kg
Weight of case	6.34 kg
EAN of inner box	4038917234252
EAN case	4038917234245

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**