



Product Specification
**Blood gas capillaries, calcium-balanced lithium
heparin, 125 µl**

Page 1



Product description

Order number	19.930.125
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	80 IU/ml blood
Application	blood gas analysis

Size

Sample volume	125 µl
Diameter	2.05 mm
Length of product	110 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Blood gas capillaries, calcium-balanced lithium
heparin, 125 µl**

Page 2

Material & colours

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	250
Type of smallest subpackaging	container
Piece(s) / inner box	250
Piece (s) / outer case	16000
Piece(s) / pallet	1280000
Depth of case	378 mm
Width of case	300 mm
Height of case	180 mm
Case volume	0.0204 cbm
Weight of product	0.0003 kg
Weight of case	6.64 kg
EAN of inner box	4038917234283
EAN case	4038917234276

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT