



Product Specification
**Blood gas capillaries, calcium-balanced lithium
heparin, 125 µl**

Page 1



Product description

| | |
|---------------------|--|
| Order number | 19.930.125 |
| Product description | Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 125 µl, (LxØ): 110 x 2.05 mm, 250 piece(s)/container |

Product characteristics

| | |
|---------------------------|----------------------------------|
| Type of collection | capillary |
| Type of preparation | calcium-balanced lithium heparin |
| Preparation concentration | 70 IU/ml blood |
| Application | blood gas analysis |

Size

| | |
|-------------------|---------|
| Sample volume | 125 µl |
| Diameter | 2.05 mm |
| Length of product | 110 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polyethylene therephthalate (PET) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Latex-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 250 |
| Type of smallest subpackaging | container |
| Piece(s) / inner box | 250 |
| Piece (s) / outer case | 16000 |
| Piece(s) / pallet | 1280000 |
| Depth of case | 378 mm |
| Width of case | 288 mm |
| Height of case | 180 mm |
| Case volume | 0.0196 cbm |
| Weight of product | 0.0003 kg |
| Weight of case | 6.64 kg |
| EAN of inner box | 4038917234283 |
| EAN case | 4038917234276 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

