



Product Specification
**Blood gas capillaries, calcium-balanced lithium
heparin, 140 µl**

Page 1



Product description

Order number	19.930.140
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 140 µl, (LxØ): 75 x 2.3 mm, 200 piece(s)/container

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	80 IU/ml blood
Application	blood gas analysis

Size

Sample volume	140 µl
Diameter	2.3 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Blood gas capillaries, calcium-balanced lithium
heparin, 140 µl**

Page 2

Material & colours

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	container
Piece(s) / inner box	200
Piece (s) / outer case	12800
Piece(s) / pallet	921600
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0002 kg
Weight of case	4.22 kg
EAN of inner box	4038917234313
EAN case	4038917234306

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT