



## Product description

Order number	51.931.100
Product description	Blood gas capillary set, preparation: calcium-balanced lithium heparin, nominal volume: 100 µl, (LxØ): 100 x 2.05 mm, set components: 1,000 PET plastic capillaries, 100 µl (4 x 250 pieces REF 19.930.100), 2,000 white quick-release caps (4 x 500 pieces REF 65.935.205) and 1,000 mixing wires (4 x 250 pieces REF 95.936), 1 piece(s)/case

## Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	80 IU/ml blood
Application	blood gas analysis

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Blood gas capillary set, calcium-balanced lithium  
heparin, 100 µl**

**Page 2**

**Size**

Sample volume	100 µl
Diameter	2.05 mm
Length of product	100 mm

**Material & colours**

Product material	Polyethylene therephthalate (PET)
Closure material	Thermoplastic elastomers (TPE)

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	432
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.816 kg
Weight of case	17.11 kg
EAN of inner box	4038917234375
EAN case	4038917234368

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**