



Product description

Order number	51.931.140
Product description	Blood gas capillary set, preparation: calcium-balanced lithium heparin, nominal volume: 140 µl, (LxØ): 75 x 2.3 mm, set components: 800 PET plastic capillaries, 140 µl (4 x 200 pieces REF 19.930.140), 1,600 blue quick-release caps (4 x 400 pieces REF 65.935.230) and 1,000 mixing wires (4 x 250 pieces REF 95.936), 1 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	80 IU/ml blood
Application	blood gas analysis

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Blood gas capillary set, calcium-balanced lithium
heparin, 140 µl**

Page 2

Size

Sample volume	140 µl
Diameter	2.3 mm
Length of product	75 mm

Material & colours

Product material	Polyethylene therephthalate (PET)
Closure material	Thermoplastic elastomers (TPE)

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	432
Depth of case	778 mm
Width of case	368 mm
Height of case	195 mm
Case volume	0.0558 cbm
Weight of product	0.5947 kg
Weight of case	13.24 kg
EAN of inner box	4038917234436
EAN case	4038917234429

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

