



**Product Specification**  
**Quick-release cap 1.60?2.10 mm for blood gas capillaries 100 µl / 125 µl, white**

**Page 1**



**Product description**

Order number	65.935.205
Product description	Quick-release cap 1.60?2.10 mm for blood gas capillaries 100 µl / 125 µl, white, blood gas analysis, 500 piece(s)/bag

**Product characteristics**

Type of collection	capillary
Application	blood gas analysis
Suitable for	blood gas capillaries 100 µl / 125 µl

**Material & colours**

Product material	Thermoplastic elastomers (TPE)
Colour of product	white

**Purity & certification**

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Quick-release cap 1.60?2.10 mm for blood gas**  
**capillaries 100 µl / 125 µl, white**

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	32000
Piece(s) / pallet	2560000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.0002 kg
Weight of case	8.42 kg
EAN of inner box	4038917234559
EAN case	4038917234542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

