



Product description

Order number	65.935.230
Product description	Quick-release cap 1.60?2.10 mm for blood gas capillaries 140 µl / 175 µl, blue, blood gas analysis, 400 piece(s)/bag

Product characteristics

Type of collection	capillary
Application	blood gas analysis
Suitable for	blood gas capillaries 140 µl / 175 µl

Material & colours

Product material	Thermoplastic elastomers (TPE)
Colour of product	blue

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Quick-release cap 1.60?2.10 mm for blood gas
capillaries 140 µl / 175 µl, blue

Page 2

Packaging

Minimum order qty.	400
Type of smallest subpackaging	bag
Piece(s) / inner box	400
Piece (s) / outer case	32000
Piece(s) / pallet	2048000
Depth of case	398 mm
Width of case	298 mm
Height of case	215 mm
Case volume	0.0255 cbm
Weight of product	0.0002 kg
Weight of case	8.13 kg
EAN of inner box	4038917234580
EAN case	4038917234573

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT