



**Product Specification**  
**Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm**

**Page 1**



**Product description**

Order number	95.937
Product description	Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm, blood gas analysis, 10 piece(s)/case

**Product characteristics**

Type of collection	capillary
Application	blood gas analysis
Suitable for	blood gas capillaries

**Size**

Diameter	19 mm
Height of product	25 mm
Length of product	25 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Magnet for blood gas capillaries, red, (LxØ): 25 x 19 mm**

**Page 2**

**Material & colours**

Colour of product red

**Purity & certification**

Product category not regulated

CE certified no CE mark

Batched yes

**Packaging**

Minimum order qty. 10

Type of smallest subpackaging case

Piece(s) / inner box 10

Piece (s) / outer case 160

Piece(s) / pallet 51840

Depth of case 268 mm

Width of case 191 mm

Height of case 100 mm

Case volume 0.0051 cbm

Weight of product 0.0072 kg

Weight of case 1.29 kg

EAN of inner box 4038917234641

EAN case 4038917234634

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**