



Product Specification
S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with paper label

Page 1



Product description

Order number	05.1079.001
Product description	S-Monovette® , blood sedimentation, 4NC, preparation: Trisodium citrate solution, nominal volume: 2 ml, (LxØ) with cap: 82 x 11 mm, with paper label, label / print: violet, cap violet, colour code EU, sterile, 50 piece(s)/case

Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate solution
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with paper label
Colour of print/label	violet
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Application	blood sedimentation
Suitable for	Sediplus® S 100

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with paper label

Page 2

Size

Sample volume	2 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

Material & colours

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with paper label

Page 3

Case volume	0,0126 cbm
Weight of product	0,0041 kg
Weight of case	2,58 kg
EAN of inner box	4038917125499
EAN case	4038917020084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT