



Product Specification
S-Monovette® , ESR, 2 ml, cap black, (LxØ): 66 x 11 mm, with paper label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 05.1079.100 |
| Product description | S-Monovette®, blood sedimentation, 4NC, preparation: Trisodium citrate solution, nominal volume: 2 ml, (LxØ) with cap: 82 x 11 mm, with paper label, label / print: black, cap black, colour code ISO, sterile, 50 piece(s)/case |

Product characteristics

| | |
|---------------------------|----------------------------|
| Type of collection | venous |
| Type of preparation | Trisodium citrate solution |
| Preparation concentration | 0.106 mol/l |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | black |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Type of cap | screw cap |
| Application | blood sedimentation |
| Suitable for | Sediplus® S 100 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® , ESR, 2 ml, cap black, (LxØ): 66 x 11 mm, with paper label

Size

| | |
|----------------------------------|--------|
| Sample volume | 2 ml |
| Diameter | 11 mm |
| Length including cap | 82 mm |
| Length excluding cap | 66 mm |
| Length including cap and plunger | 103 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Cap material | High Density Polyethylene (HD-PE) |
| Colour of cap | black |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In-vitro diagnostic, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------|--------|
| Minimum order qty. | 500 |
| Type of inner packaging | case |
| Piece(s) / inner box | 50 |
| Piece(s) / case | 500 |
| Piece(s) / pallet | 56000 |
| Depth of box | 145 mm |
| Width of box | 108 mm |
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
S-Monovette® , ESR, 2 ml, cap black, (LxØ): 66 x 11 mm, with paper label

| | |
|-------------------|---------------|
| Case volume | 0,0126 cbm |
| Weight of product | 0,0041 kg |
| Weight of case | 2,58 kg |
| EAN of inner box | 4038917125505 |
| EAN case | 4038917020114 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

