



**Product Specification**  
**S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with plastic label**

**Page 1**



### Product description

Order number	05.1079
Product description	S-Monovette® , blood sedimentation, 4NC, preparation: Trisodium citrate solution, nominal volume: 2 ml, (LxØ) with cap: 82 x 11 mm, with plastic label, label / print: transparent/violet, cap violet, colour code EU, sterile, 50 piece(s)/case

### Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate solution
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/violet
Graduation	fill mark
Cap	membrane screw cap
Type of cap	screw cap
Application	blood sedimentation
Suitable for	Sediplus® S 100

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-01-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with plastic label**

**Size**

Sample volume	2 ml
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

**Material & colours**

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with plastic label**

Case volume	0,0126 cbm
Weight of product	0,004 kg
Weight of case	2,56 kg
EAN of inner box	4038917085298
EAN case	4038917003995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

