



Product description

| | |
|---------------------|--|
| Order number | 05.1079 |
| Product description | S-Monovette® , blood sedimentation, 4NC, preparation: Citrate, 4NC, nominal volume: 2 ml, (LxØ) with cap: 82 x 11 mm, with plastic label, label/print: transparent/violet, cap violet, colour code EU, sterile, 50 piece(s)/case |

Product characteristics

| | |
|---------------------------|---------------------|
| Type of collection | venous |
| Type of preparation | Citrate, 4NC |
| Preparation concentration | 0.106 mol/l |
| Colour code | EU |
| Label/ Print | with plastic label |
| Colour of print/label | transparent/violet |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Application | blood sedimentation |
| Suitable for | Sediplus® S 100 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|----------------------------------|--------|
| Sample volume | 2 ml |
| Diameter | 11 mm |
| Length including cap | 82 mm |
| Length excluding cap | 66 mm |
| Length including cap and plunger | 103 mm |

Material & colours

| | |
|-------------------------|-----------------------------------|
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | violet |
| Membrane material | rubber |
| Tube material | Polypropylene (PP) |
| Piston material | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS) |

Purity & certification

| | |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR |
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Purity standard | sterile |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|-------------------------------|--------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 56000 |
| Depth of box | 145 mm |
| Width of box | 108 mm |
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® , ESR, 2 ml, cap violet, (LxØ): 66 x 11 mm, with plastic label

| | |
|-------------------|---------------|
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.004 kg |
| Weight of case | 2.56 kg |
| EAN of inner box | 4038917085298 |
| EAN case | 4038917003995 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

