



**Product Specification**  
**S-Sedivette® , 3.5 ml, cap violet, (LxØ): 130 x 8 mm, with paper label**

**Page 1**



### Product description

Order number	06.1690.001
Product description	S-Sedivette® , blood sedimentation, 4NC, preparation: Trisodium citrate/citric acid buffer solution, nominal volume: 3.5 ml, (LxØ) with cap: 146 x 8 mm, with paper label, label/print: violet, cap violet, colour code EU, sterile, 50 piece(s)/case

### Product characteristics

Type of collection	venous
Type of preparation	Trisodium citrate/citric acid buffer solution
Preparation concentration	0.105 mol/l pH 5.5
Colour code	EU
Label/ Print	with paper label
Colour of print/label	violet
Graduation	no
Cap	membrane screw cap
Type of cap	screw cap
Application	blood sedimentation
Suitable for	Sediplus® S 2000 NX

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Sedivette® , 3.5 ml, cap violet, (LxØ): 130 x 8 mm, with paper label**

**Page 2**

**Size**

Sample volume	3.5 ml
Diameter	8 mm
Length including cap	146 mm
Length excluding cap	130 mm
Length including cap and plunger	167 mm

**Material & colours**

Colour of product	transparent
Cap material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of inner packaging	case
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	36000
Depth of box	143 mm
Width of box	171 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Sedivette® , 3.5 ml, cap violet, (LxØ): 130 x 8 mm, with paper label**

**Page 3**

Case volume	0.0192 cbm
Weight of product	0.0052 kg
Weight of case	3.29 kg
EAN of inner box	4038917133111
EAN case	4038917024709

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**