



### Product description

Order number	18.1325
Product description	Microvette® CB 200 NC, capillary blood collection, preparation: ESR/Citrate Buffer (1:5), nominal volume: 200 µl, (LxØ) with cap: 46 x 10.8 mm, cap violet, push cap, flat base, incl. sedimentation pipette, 50 piece(s)/bag

### Product characteristics

Type of collection	capillary
Type of preparation	ESR/Citrate Buffer (1:5)
Colour code	EU
Type of cap	push cap
Application	capillary blood collection
Base shape	flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





## Size

Sample volume	200 µl
Diameter	10.8 mm
Length of product	40 mm
Length including cap	46 mm
Length excluding cap	40 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Low Density Polyethylene (LD-PE)
Colour of cap	violet

## Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	32000
Depth of box	271 mm
Width of box	126 mm
Height of box	36 mm
Depth of case	283 mm
Width of case	268 mm
Height of case	210 mm
Case volume	0.0159 cbm
Weight of product	0.0036 kg
Weight of case	2.6 kg
EAN of inner box	4038917086202
EAN case	4038917004787

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

