



Product description

Order number 86.1996
Product description ESR pipette, with fitted piston insert

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1000
Type of inner packaging case
Piece(s) / inner box 100
Piece(s) / case 1000
Piece(s) / pallet 42000
Depth of box 233 mm
Width of box 164 mm
Height of box 76 mm
Depth of case 388 mm
Width of case 338 mm
Height of case 250 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR pipette

Page 2

Case volume	0,0328 cbm
Weight of product	0,0076 kg
Weight of case	8,61 kg
EAN of inner box	4038917088138
EAN case	4038917006231

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT