



Product description

| | |
|---------------------|--|
| Order number | 86.1996 |
| Product description | ESR pipette, with fitted piston insert |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | In-vitro diagnostic, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------|--------|
| Minimum order qty. | 1000 |
| Type of inner packaging | case |
| Piece(s) / inner box | 100 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 42000 |
| Depth of box | 233 mm |
| Width of box | 164 mm |
| Height of box | 76 mm |
| Depth of case | 388 mm |
| Width of case | 338 mm |
| Height of case | 250 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR pipette

Page 2

| | |
|-------------------|---------------|
| Case volume | 0.0328 cbm |
| Weight of product | 0.0076 kg |
| Weight of case | 8.61 kg |
| EAN of inner box | 4038917088138 |
| EAN case | 4038917006231 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT