



Product description

Order number 90.1060.062
Product description ESR rack, suitable for S-Monovette® , without back panel

Product characteristics

Design without back panel
Suitable for S-Monovette®

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 20
Piece(s) / pallet 1560
Depth of case 428 mm
Width of case 323 mm
Height of case 130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR rack, suitable for S-Monovette®

Page 2

Case volume	0.018 cbm
Weight of product	1.167 kg
Weight of case	23.8 kg
EAN of inner box	4038917124935
EAN case	4038917019545

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-10-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT