



Product description

Order number	90.1060
Product description	ESR rack, suitable for S-Monovette® , with back panel

Product characteristics

Design	with back panel
Suitable for	S-Monovette®

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece (s) / outer case	8
Piece(s) / pallet	160
Depth of case	598 mm
Width of case	400 mm
Height of case	350 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
ESR rack, suitable for S-Monovette®

Page 2

Case volume	0.0837 cbm
Weight of product	0.5005 kg
Weight of case	7.1 kg
EAN of inner box	4038917085205
EAN case	4038917003896

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-08

This document was prepared by EDP support and is valid without signature.



SARSTEDT