



Product description

Order number 90.1060
Product description ESR rack, suitable for S-Monovette® , with back panel

Product characteristics

Design with back panel
Suitable for S-Monovette®

Purity & certification

Product category In vitro diagnostic device, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece (s) / outer case 8
Piece(s) / pallet 160
Depth of case 263 mm
Width of case 98 mm
Height of case 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
ESR rack, suitable for S-Monovette®

Page 2

Case volume	0.0086 cbm
Weight of product	0.5005 kg
Weight of case	7.1 kg
EAN of inner box	4038917085205
EAN case	4038917003896

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT