



## Product description

Order number 90.1060  
Product description ESR rack, suitable for S-Monovette® , with back panel

## Product characteristics

Design with back panel  
Suitable for S-Monovette®

## Purity & certification

Product category In-vitro diagnostic, CE  
CE certified CE  
Batched yes

## Packaging

Minimum order qty. 1  
Piece(s) / inner box 1  
Piece(s) / case 8  
Piece(s) / pallet 160  
Depth of case 263 mm  
Width of case 98 mm  
Height of case 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**ESR rack, suitable for S-Monovette®**

**Page 2**

Case volume	0.0086 cbm
Weight of product	0.5005 kg
Weight of case	7.1 kg
EAN of inner box	4038917085205
EAN case	4038917003896

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-10-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**