



Product Specification
ESR rack, suitable for S-Sedivette® , 10 measuring stations

Page 1



Product description

Order number 90.1090
Product description ESR rack, suitable for S-Sedivette® , 10 measuring stations

Product characteristics

Suitable for S-Sedivette®
Number of positions 10

Material & colours

Product material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	320
Depth of box	73 mm
Width of box	73 mm
Height of box	211 mm
Depth of case	598 mm
Width of case	400 mm
Height of case	350 mm
Case volume	0.0837 cbm
Weight of product	0.6018 kg
Weight of case	11.75 kg
EAN of inner box	4038917085304
EAN case	4038917004008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

