



Product Specification
ESR rack, suitable for S-Sedivette® , 10 measuring stations

Page 1



Product description

Order number	90.1090
Product description	ESR rack, suitable for S-Sedivette® , 10 measuring stations

Product characteristics

Suitable for	S-Sedivette®
Number of positions	10

Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	16
Piece(s) / pallet	320
Depth of case	598 mm
Width of case	398 mm
Height of case	350 mm
Case volume	0.0833 cbm
Weight of product	0.6018 kg
Weight of case	11.75 kg
EAN of inner box	4038917085304
EAN case	4038917004008

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

