




Product Specification
ESR rack, suitable for Microvette® , 10 measuring stations

Page 1



Product description

Order number 90.1091
Product description ESR rack, suitable for Microvette® , 10 measuring stations

Product characteristics

Suitable for Microvette®
Number of positions 10

Material & colours

Product material Polypropylene (PP)

Purity & certification

Product category In-vitro diagnostic (IVD)
CE certified CE - manufacturer's self-declaration
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	8
Piece(s) / pallet	160
Depth of case	598 mm
Width of case	400 mm
Height of case	350 mm
Case volume	0.0837 cbm
Weight of product	0.2505 kg
Weight of case	5.1 kg
EAN of inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

