



Product description

Order number 90.1091
Product description ESR rack, suitable for Microvette®

Product characteristics

Suitable for Microvette®

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Batched yes

Packaging

Minimum order qty. 1
Type of inner packaging case
Piece(s) / inner box 1
Piece(s) / case 8
Piece(s) / pallet 160
Depth of case 263 mm
Width of case 98 mm
Height of case 335 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
ESR rack, suitable for Microvette®

Page 2

Case volume	0.0086 cbm
Weight of product	0.2505 kg
Weight of case	5.1 kg
EAN of inner box	4038917085311
EAN case	4038917004015

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-07-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT