




Product Specification
ESR rack, suitable for Microvette® , 10 measuring stations

Page 1



Product description

Order number 90.1091
Product description ESR rack, suitable for Microvette® , 10 measuring stations

Product characteristics

Suitable for Microvette®
Number of positions 10

Material & colours

Product material Polypropylene (PP)

Purity & certification

Product category In vitro diagnostic device, CE
CE certified CE
Batched yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 8 |
| Piece(s) / pallet | 160 |
| Depth of case | 263 mm |
| Width of case | 98 mm |
| Height of case | 335 mm |
| Case volume | 0.0086 cbm |
| Weight of product | 0.2505 kg |
| Weight of case | 5.1 kg |
| EAN of inner box | 4038917085311 |
| EAN case | 4038917004015 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

