



Product description

Order number 90.1092.730
Product description Barcode scanner, for Sediplus® S 2000 NX

Product characteristics

Suitable for Sediplus® S 2000 NX

Size

Width of product 85 mm
Height of product 60 mm
Length of product 150 mm

Purity & certification

Product category no medical device | IVD
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	1
Piece(s) / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.366 kg
Weight of case	0.47 kg
EAN of inner box	4038917927130
EAN case	4038917927147

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

