



Product description

Order number	90.188.055
Product description	Paper roll, for thermo printer item no 90.189.720 & 90.1092.720

Product characteristics

Suitable for	thermo printer item no 90.189.720 & 90.1092.720
--------------	---

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Case volume	0.0007 cbm
Weight of product	0.404 kg
Weight of case	0.43 kg
EAN of inner box	4038917102759

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
**Paper roll, for thermo printer item no 90.189.720 &
90.1092.720**

EAN case

4038917008884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT