



Product Specification
Sediplus® , automatic 40-channel blood sedimentation measuring device, 40 measuring stations



Product description

Order number 90.189.700
Product description Sediplus® , S 2000, automatic 40-channel blood sedimentation measuring device, 40 measuring stations

Product characteristics

Number of positions 40
Supply voltage 230 V

Purity & certification

Product category IVD product
CE certified CE
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Sediplus® , automatic 40-channel blood sedimentation measuring device, 40 measuring stations

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	1
Case volume	0,0371 cbm
Weight of product	6,303 kg
Weight of case	7,6 kg
EAN of inner box	4038917102957
EAN case	4038917009003

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

