



## Product Specification

**Sediplus® , automatic 40-channel blood sedimentation measuring device, 40 measuring stations, 230 V**

**Page 1**



### Product description

Order number	90.189.700
Product description	Sediplus® , S 2000, automatic 40-channel blood sedimentation measuring device, 40 measuring stations, supply voltage: 230 V

### Product characteristics

Number of positions	40
Supply voltage	230 V

### Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Batched	no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-09-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sediplus® , automatic 40-channel blood sedimentation measuring device, 40 measuring stations, 230 V**

**Page 2**

**Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	1
Case volume	0.0371 cbm
Weight of product	6.303 kg
Weight of case	7.6 kg
EAN of inner box	4038917102957
EAN case	4038917009003

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-09-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**