



### **Product description**

Order number 90.189.730  
Product description Barcode scanner, for Sediplus® S100/S200/S2000

### **Product characteristics**

Suitable for Sediplus® S100/S200/S2000

### **Purity & certification**

Product category no medical device | IVD  
CE certified CE  
Batched no

### **Packaging**

Minimum order qty. 1  
Piece(s) / inner box 1  
Piece(s) / case 1  
Case volume 0,0026 cbm  
Weight of product 0,35 kg  
Weight of case 0,4 kg  
EAN of inner box 4038917103039  
EAN case 4038917009034

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

