




Product Specification

Sediplus® , 10-channel blood sedimentation measuring device, incl. mains adapter 12 V DC, 10 measuring stations

Page 1



Product description

Order number 90.189.800
Product description Sediplus® , S 100, 10-channel blood sedimentation measuring device, incl. mains adapter 12 V DC, 10 measuring stations

Product characteristics

Number of positions 10
Supply voltage 12 V

Purity & certification

Product category In-vitro diagnostic, CE
CE certified CE
Batched no

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Sediplus® , 10-channel blood sedimentation
measuring device, incl. mains adapter 12 V DC, 10
measuring stations**

Page 2

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	36
Case volume	0,0342 cbm
Weight of product	4,8995 kg
Weight of case	4,9 kg
EAN of inner box	4038917226400
EAN case	4038917226394

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT