



Product description

Order number	85.1015
Product description	Safety lancet, Mini, Ø needle: 28 G, light blue, penetration depth: 1.6 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

Product characteristics

Application	capillary blood collection
Example of use	Young children, adults: diabetic patients, oncology, determining blood sugar
Puncture site	finger, ear, heel
Penetration depth	1.6 mm
Safety product	yes

Size

Blood volume	0,1 µl - 10 µl
Diameter of needle	28 G
Width of product	21 mm
Length of product	52 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Safety lancet, Mini, Ø needle: 28 G, penetration
depth: 1.6 mm**

Page 2

Material & colours

Colour of product light blue

Purity & certification

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order qty.	1200
Type of smallest subpackaging	case
Piece(s) / inner box	200
Piece (s) / outer case	1200
Piece(s) / pallet	108000
Depth of box	118 mm
Width of box	118 mm
Height of box	151 mm
Depth of case	366 mm
Width of case	248 mm
Height of case	173 mm
Case volume	0.0157 cbm
Weight of product	0.0028 kg
Weight of case	4.08 kg
EAN of inner box	4038917085021
EAN case	4038917057622

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-12-14

This document was prepared by EDP support and is valid without signature.



SARSTEDT