



### Product description

Order number	85.1016
Product description	Safety lancet, Normal, Ø needle: 21 G, green, penetration depth: 1.8 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

### Product characteristics

Application	capillary blood collection
Example of use	Children, adults: haematocrit, standard blood count
Puncture site	finger, ear, heel
Penetration depth	1.8 mm
Safety product	yes

### Size

Blood volume	10 µl - 50 µl
Diameter of needle	21 G
Width of product	21 mm
Length of product	52 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Colour of product green

### **Purity & certification**

Product category	Medical device (MD)
CE certified	CE0124 - DEKRA
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1200
Type of smallest subpackaging	case
Piece(s) / inner box	200
Piece (s) / outer case	1200
Piece(s) / pallet	108000
Depth of box	118 mm
Width of box	118 mm
Height of box	151 mm
Depth of case	366 mm
Width of case	248 mm
Height of case	173 mm
Case volume	0.0157 cbm
Weight of product	0.0028 kg
Weight of case	4.08 kg
EAN of inner box	4038917085038
EAN case	4038917057639

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

