



### Product description

Order number	85.1017
Product description	Safety lancet, Extra, Ø needle: 18 G, yellow, penetration depth: 1.8 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

### Product characteristics

Application	capillary blood collection
Example of use	Children, Adults, obese adults: capillary blood gas analysis, CBC
Puncture site	finger, ear, heel
Penetration depth	1.8 mm
Safety product	yes

### Size

Blood volume	20 µl - 100 µl
Diameter of needle	18 G
Width of product	21 mm
Length of product	52 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Safety lancet, Extra, Ø needle: 18 G, penetration  
depth: 1.8 mm**

**Page 2**

**Material & colours**

Colour of product                      yellow

**Purity & certification**

Product category                      Medical device  
CE certified                              CE 0124 (DEKRA)  
Sterilisation                              Gamma irradiation  
Batched                                    yes

**Packaging**

Minimum order qty.                    1200  
Type of inner packaging               case  
Piece(s) / inner box                   200  
Piece(s) / case                         1200  
Piece(s) / pallet                        108000  
Depth of box                            118 mm  
Width of box                             118 mm  
Height of box                            151 mm  
Depth of case                            366 mm  
Width of case                            248 mm  
Height of case                           173 mm  
Case volume                             0.0157 cbm  
Weight of product                      0.0029 kg  
Weight of case                           4.08 kg  
EAN of inner box                        4038917085045  
EAN case                                 4038917057646

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2022-05-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**