



Product description

Order number	85.1018
Product description	Safety lancet, Super, violet, penetration depth: 1.6 mm, blade width: 1.5 mm, for capillary blood collection, for puncturing finger, ear, heel, 200 piece(s)/case

Product characteristics

Application	capillary blood collection
Example of use	Children, Adults when very large blood volumes are required, neonates: CBC, reticulocytes, metabolism test
Puncture site	finger, ear, heel
Penetration depth	1.6 mm
Safety product	yes

Size

Blood volume	100 µl - 500 µl
Width of product	21 mm
Length of product	52 mm
Width of blade	1.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product violet

Purity & certification

Product category Medical device
CE certified CE 0124 (DEKRA)
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1200
Type of inner packaging case
Piece(s) / inner box 200
Piece(s) / case 1200
Piece(s) / pallet 108000
Depth of box 118 mm
Width of box 118 mm
Height of box 151 mm
Depth of case 366 mm
Width of case 248 mm
Height of case 173 mm
Case volume 0.0157 cbm
Weight of product 0.0029 kg
Weight of case 4.08 kg
EAN of inner box 4038917085052
EAN case 4038917057653

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

