



Product description

Order number	85.1019
Product description	Safety lancet, Neonatal, pink, penetration depth: 1.2 mm, blade width: 1.5 mm, for capillary blood collection, for puncturing heel, finger, 200 piece(s)/case

Product characteristics

Application	capillary blood collection
Example of use	Children, Adults when large blood volumes are required, neonates: standard blood count, blood sugar, bilirubin
Puncture site	heel, finger
Penetration depth	1.2 mm
Safety product	yes

Size

Blood volume	50 µl - 250 µl
Width of product	21 mm
Length of product	52 mm
Width of blade	1.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Colour of product pink

Purity & certification

Product category Medical device (MD)
CE certified CE0124 - DEKRA
Purity standard sterile
Sterilisation Gamma irradiation
Batched yes

Packaging

Minimum order qty. 1200
Type of smallest case
subpackaging
Piece(s) / inner box 200
Piece (s) / outer case 1200
Piece(s) / pallet 108000
Depth of box 118 mm
Width of box 118 mm
Height of box 151 mm
Depth of case 366 mm
Width of case 248 mm
Height of case 173 mm
Case volume 0.0157 cbm
Weight of product 0.0029 kg
Weight of case 4.08 kg
EAN of inner box 4038917085069
EAN case 4038917057660

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

