



### **Product description**

Order number	85.1026
Product description	Incision safety lancet, Safety-Heel® Newborn, light blue/white, penetration depth: 1 mm, blade width: 2.5 mm, for heel puncture, 50 piece(s)/case

### **Product characteristics**

Example of use	Blood collection from newborn babies
Puncture site	Heel
Penetration depth	1 mm
Safety product	yes

### **Size**

Width of product	35 mm
Length of product	40 mm
Width of blade	2.5 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Incision safety lancet, Safety-Heel® Newborn,**  
**penetration depth: 1 mm**

**Page 2**

**Material & colours**

Colour of product                      light blue/white

**Purity & certification**

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

**Packaging**

Minimum order qty.	200
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	200
Piece(s) / pallet	24000
Depth of box	255 mm
Width of box	120 mm
Height of box	70 mm
Depth of case	273 mm
Width of case	262 mm
Height of case	160 mm
Case volume	0.0114 cbm
Weight of product	0.0094 kg
Weight of case	2.08 kg
EAN of inner box	6945630110474
EAN case	6945630110535

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2025-12-14

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**