



**Product Specification**  
**Multivette® 600 EDTA K3E, 600 µl, cap violet,  
screw cap**

**Page 1**



**Product description**

Order number	15.1671.100
Product description	Multivette® 600 EDTA K3E, venous and capillary blood samples, preparation: K3 EDTA, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap violet, screw cap, flat base, 100 piece(s)/case

**Product characteristics**

Type of collection	Venous and capillary
Type of preparation	K3 EDTA
Colour code	ISO
Label/ Print	with print
Colour of print/label	violet
Closure type	screw cap
Application	venous and capillary blood samples
Base shape	flat base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Multivette® 600 EDTA K3E, 600 µl, cap violet,  
screw cap**

**Page 2**

**Size**

Sample volume	600 µl
Diameter	8 mm
Length including cap	47.9 mm
Length excluding cap	44.4 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	63000
Depth of box	201 mm
Width of box	81 mm
Height of box	81 mm
Depth of case	423 mm
Width of case	213 mm
Height of case	185 mm
Case volume	0.0167 cbm
Weight of product	0.0022 kg
Weight of case	2.79 kg
EAN of inner box	4038917132954
EAN case	4038917024594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**