



**Product Specification**  
**Microvette® CB 300 Lithium heparin LH, 300 µl,  
cap orange, push cap, flat base**

**Page 1**



### Product description

Order number	16.443
Product description	Microvette® CB 300 Lithium heparin LH, Veterinary medicine, preparation: Lithium heparin, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap orange, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag

### Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin
Colour code	EU
Cap	Soft stopper with retaining loop
Closure type	push cap
Application	Veterinary medicine
Base shape	flat base
Recommended centrifugation	2000 x g
Centrifugation time	5 Min
Shape of inner tube	trumpet-shaped

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-09-07

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



### **Size**

Sample volume	300 µl
Diameter	10.8 mm
Length including cap	46.5 mm
Length excluding cap	43.7 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	orange

### **Purity & certification**

Product category	no medical device   IVD
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	80000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0.0196 cbm
Weight of product	0.002 kg
Weight of case	2.35 kg
EAN of inner box	4038917268295
EAN case	4038917282727

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

