



Product description

Order number	16.443
Product description	Microvette® CB 300 Lithium heparin, capillary blood collection, preparation: Lithium heparin, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap orange, push cap, flat base, inner receptacle: trumpet-shaped

Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin
Colour code	EU
Type of cap	push cap
Application	capillary blood collection
Base shape	flat base
Recommended centrifugation	2000 U/min
Centrifugation time	5 Min
Shape of inner tube	trumpet-shaped

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	300 µl
Diameter	10,8 mm
Length including cap	46,5 mm
Length excluding cap	43,7 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Low Density Polyethylene (LD-PE)
Colour of cap	orange

Purity & certification

Product in accordance with MDD or IVDD	IVD product
CE certified	CE
Lot or serial number	yes

Packaging

Minimum order qty.	1000
Type of inner packaging	bag
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Depth of case	378 mm
Width of case	288 mm
Height of case	180 mm
Case volume	0,0196 cbm
Weight of product	0,002 kg
Weight of case	2,35 kg
EAN of inner box	4038917268295
EAN case	4038917282727

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

