



**Product description**

Order number	17.2111.050
Product description	Minivette® POCT neutral Z, nominal volume: 50 µl, without preparation, white, 200 piece(s)/bag, 200 piece(s)/case

**Product characteristics**

Type of collection	capillary
Type of preparation	without preparation
Application	capillary blood collection
Coding colour	white

**Size**

Sample volume	50 µl
---------------	-------

**Material & colours**

Product material	Cellulose acetate (CA)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	200
Piece (s) / outer case	2000
Piece(s) / pallet	144000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.001 kg
Weight of case	2.66 kg
EAN of inner box	4038917310314
EAN case	4038917310307

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

