



Product description

| | |
|---------------------|--|
| Order number | 17.2111.200 |
| Product description | Minivette® POCT neutral Z, nominal volume: 200 µl, without preparation, white, 150 piece(s)/bag, 150 piece(s)/case |

Product characteristics

| | |
|---------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | without preparation |
| Application | capillary blood collection |
| Coding colour | white |

Size

| | |
|---------------|--------|
| Sample volume | 200 µl |
|---------------|--------|

Material & colours

| | |
|-------------------|------------------------|
| Product material | Cellulose acetate (CA) |
| Colour of product | transparent |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD) |
| CE certified | CE - manufacturer's self-declaration |
| Latex-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 150 |
| Piece (s) / outer case | 1500 |
| Piece(s) / pallet | 108000 |
| Depth of box | 140 mm |
| Width of box | 99 mm |
| Height of box | 117 mm |
| Depth of case | 513 mm |
| Width of case | 293 mm |
| Height of case | 140 mm |
| Case volume | 0.021 cbm |
| Weight of product | 0.0011 kg |
| Weight of case | 2.34 kg |
| EAN of inner box | 4038917344005 |
| EAN case | 4038917343992 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

