



Product Specification

**Minivette® POCT Lithium heparin, 200 µl, plunger
orange, colour code EU, 150 piece(s)/bag**

Page 1



Product description

Order number	17.2112.200
Product description	Minivette® POCT Lithium heparin, nominal volume: 200 µl, preparation: Heparin, orange, colour code EU, 150 piece(s)/bag, 150 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Heparin
Colour code	EU
Application	capillary blood collection
Coding colour	orange

Size

Sample volume	200 µl
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-19

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Material & colours

Product material	Cellulose acetate (CA)
Colour of product	transparent

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	1500
Type of inner packaging	bag
Piece(s) / inner box	150
Piece(s) / case	1500
Piece(s) / pallet	108000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0,021 cbm
Weight of product	0,0011 kg
Weight of case	2,37 kg
EAN of inner box	4038917337052
EAN case	4038917337045

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

