



**Product Specification**  
**End-to-end capillary, Ammonium heparin AH, 10  
µl, Ø: 1.25 mm, glass**

**Page 1**

**Product description**

Order number	19.414.001
Product description	End-to-end capillary, Ammonium heparin AH, capillary blood collection, preparation: Ammonium heparin, nominal volume: 10 µl, Ø: 1.25 mm, glass, 100 piece(s)/container

**Product characteristics**

Type of collection	capillary
Type of preparation	Ammonium heparin
Application	capillary blood collection
Volume of work	10 µl

**Size**

Sample volume	10 µl
Blood volume	10 µl
Diameter	1.25 mm
Length of product	29 mm

**Material & colours**

Product material	glass
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE supplier
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	container
Piece(s) / inner box	1000
Piece (s) / outer case	50000
Piece(s) / pallet	2400000
Depth of box	117 mm
Width of box	50 mm
Height of box	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**End-to-end capillary, Ammonium heparin AH, 10**  
**µl, Ø: 1.25 mm, glass**

**Page 2**

Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0001 kg
Weight of case	8.65 kg
EAN of inner box	5706860179225
EAN case	5706860179232

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**