



Product description

| | |
|---------------------|---|
| Order number | 20.1282 |
| Product description | Microvette® 100 Lithium heparin LH, capillary blood collection, preparation: Lithium heparin, nominal volume: 100 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap orange, flat base, inner receptacle: conical, 100 piece(s)/case |

Product characteristics

| | |
|----------------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | Lithium heparin |
| Colour code | EU |
| Label/ Print | with print |
| Colour of print/label | orange |
| Cap | cap |
| Application | capillary blood collection |
| Base shape | flat base |
| Recommended centrifugation | 2000 x g |
| Centrifugation time | 5 Min |
| Shape of inner tube | conical |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Microvette® 100 Lithium heparin LH, 100 µl, cap orange, flat base

Page 2

Size

| | |
|----------------------|---------|
| Sample volume | 100 µl |
| Diameter | 10.8 mm |
| Length including cap | 47.6 mm |
| Length excluding cap | 46.6 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | orange |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 121000 |
| Depth of box | 139 mm |
| Width of box | 108 mm |
| Height of box | 66 mm |
| Depth of case | 338 mm |
| Width of case | 228 mm |
| Height of case | 160 mm |
| Case volume | 0.0123 cbm |
| Weight of product | 0.0026 kg |
| Weight of case | 3.11 kg |
| EAN of inner box | 4038917086004 |
| EAN case | 4038917057134 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

