



## Product description

|                     |  |
|---------------------|--|
| Order number        | 20.1292  |
| Product description | Microvette® 200 Lithium heparin LH, capillary blood collection, preparation: Lithium heparin LH, nominal volume: 200 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap orange, flat base, inner receptacle: conical, 100 piece(s)/case |

## Product characteristics

|                            |                            |
|----------------------------|----------------------------|
| Type of collection         | capillary                  |
| Type of preparation        | Lithium heparin LH         |
| Colour code                | EU                         |
| Label/ Print               | with print                 |
| Colour of print/label      | orange                     |
| Cap                        | cap                        |
| Application                | capillary blood collection |
| Base shape                 | flat base                  |
| Recommended centrifugation | 2000 x g                   |
| Centrifugation time        | 5 Min                      |
| Shape of inner tube        | conical                    |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

|                      |         |
|----------------------|---------|
| Sample volume        | 200 µl  |
| Diameter             | 10.8 mm |
| Length including cap | 47.6 mm |
| Length excluding cap | 46.6 mm |

## Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Product material  | Polypropylene (PP)                |
| Colour of product | transparent                       |
| Closure material  | High Density Polyethylene (HD-PE) |
| Colour of cap     | orange                            |

## Purity & certification

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

## Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 121000        |
| Depth of box                  | 139 mm        |
| Width of box                  | 108 mm        |
| Height of box                 | 66 mm         |
| Depth of case                 | 338 mm        |
| Width of case                 | 228 mm        |
| Height of case                | 160 mm        |
| Case volume                   | 0.0123 cbm    |
| Weight of product             | 0.0027 kg     |
| Weight of case                | 3.23 kg       |
| EAN of inner box              | 4038917086066 |
| EAN case                      | 4038917070386 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

