



## Product description

Order number	20.1341.102
Product description	Microvette® 500 EDTA K3E, low lead, capillary blood collection, preparation: K3 EDTA, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap violet, flat base, inner receptacle: cylindrical, 100 piece(s)/case

## Product characteristics

Type of collection	capillary
Type of preparation	K3 EDTA
Colour code	ISO
Label/ Print	with print
Colour of print/label	violet
Cap	quadruple-threaded cap
Application	capillary blood collection
Base shape	flat base
Volume of work	325 µl - 500 µl
Shape of inner tube	cylindrical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Microvette® 500 EDTA K3E, low lead, 500 µl, cap violet, flat base**

**Page 2**

**Size**

Sample volume	500 µl
Diameter	10.8 mm
Length including cap	47.6 mm
Length excluding cap	46.6 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	242000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	121000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	160 mm
Case volume	0.0123 cbm
Weight of product	0.002 kg
Weight of case	2.47 kg
EAN of inner box	4038917355377
EAN case	4038917355360

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**